



Dimension throat	
After etching	50 μm
Mask width	10 μm
Depth	20 μm
Dimension pore	
After etching diameter	90 μm
Mask width diameter	50 μm
Depth	20 μm

Note:  
The uniform network consists of 10151 pores, which are connected by throats. In this number is included the pores connected to the inlet/outlet channels



ITEM NUMBER: 02978  
DRAWN BY: WVDC  
DRAWN DATE: 22 Dec 2021  
PROJECT NO:  
EXT. REF.:

TITLE: Chip EOR UN 20.2

PROPRIETARY AND CONFIDENTIAL  
THE INFORMATION CONTAINED IN THIS DRAWING IS THE SOLE PROPERTY OF MICRONIT BV. ANY REPRODUCTION IN PART OR AS A WHOLE WITHOUT THE WRITTEN PERMISSION OF MICRONIT BV IS PROHIBITED.

UNLESS OTHERWISE SPECIFIED:  
- DIMENSIONS ARE IN MILLIMETERS  
- LINEAR TOLERANCES: ± 0.1  
- ANGULAR TOLERANCES: ± 5°  
DO NOT SCALE DRAWING

Zone	Rev	REVISION NOTE	Mod By	App By	App Date

DWG NO. 42005552  
REVISION 0  
SCALE: 2:1  
SHEET FORMAT: A3  
SHEET 1 OF 1