




DETAIL B
SCALE 50 : 1



SECTION A-A
SCALE 5 : 1



SECTION C-C
SCALE 50 : 1



DETAIL D
SCALE 50 : 1

The structures are randomly placed, resulting in throats and channels of random dimensions. The design resembles structures from cut sandstone.

The vertical crack is added to the physical rock structure, therefore around the crack, it will interact with the rock structure

Zone	1.0	First Release	IR	05 Jun 2023
Rev		REVISION NOTE	Mod By	App By

	ITEM NUMBER:	03200	DRAWN BY:	WVDC	TITLE:	Chip EOR PR 20.2 Crack 150 μm	
	THIRD ANGLE PROJECTION		DRAWN DATE:	11 May 2023	DWG NO.	42008030	REVISION
PROPRIETARY AND CONFIDENTIAL	UNLESS OTHERWISE SPECIFIED:	- DIMENSIONS ARE IN MILLIMETERS	- LINEAR TOLERANCES: ± 0.1	- ANGULAR TOLERANCES: ± 5°	MATERIAL:	Borofloat 33	1.0
DO NOT SCALE DRAWING					SCALE:	5:1	SHEET FORMAT: A3 SHEET 1 OF 1